

ARTAX

- **User Manual**

This page intentionally left blank.

Table of Contents

1	Introduction	5
1.1	Installation	5
1.2	Program Start.....	5
2	User Interface	7
2.1	Menu Structure.....	7
2.1.1	File Menu	8
2.1.2	Device Menu	8
2.1.3	Measurement Menu	9
2.1.4	Analyze Menu	9
2.1.5	Spectrum Menu	10
2.1.6	Project Menu	12
2.1.7	Options Menu	13
2.1.8	Export Menu	14
2.1.9	User Menu	15
2.1.10	Help Menu	15
2.2	Toolbar Structure.....	16
3	Working with Projects	18
3.1	Opening a Project	18
3.2	Tabs.....	19
3.3	Creating a Project	22
3.4	Copying and Saving	23
4	Working with Spectra	23
4.1	Active Spectrum.....	24
4.2	Display Options	25
4.3	Manipulating a Spectrum	29
4.4	Status Bar	30
5	Regions of Interest (ROIs)	31
6	Spectra Parameters	35
7	Exporting	36
8	Working with Methods.....	39
8.1	Method Editor	39
8.2	Measurement Panel	41
8.3	Corrections Panel.....	42
8.4	Identification Panel.....	43
8.5	Deconvolution Panel.....	44
8.6	Qualification Panel.....	44
8.7	PDZ Options Panel	44
9	Spectra Processing.....	45
9.1	Adding and Subtracting Spectra	45
9.2	Scaling a Spectrum	46
9.3	Comparing Spectra	47
10	Qualitative Analysis.....	49
10.1	Display Options	49
10.2	Deconvolution.....	51
	Appendix A: Registry Database.....	52

This page intentionally left blank.



1 Introduction

- Description** Artax is a program developed to work with S1 TITAN, TRACER 5, and CTX instruments to –
- Control system components.
 - Acquire measurement data.
 - Analyze measurement data.
- Analyzing data** Measured data can be analyzed by –
- Spectra correction (Escape, Shelf, Background).
 - Qualitative analysis.
 - Calculation of the spectral line intensities by deconvolution.
- Contacting Bruker** *Email:* support.hmp@bruker.com
Phone: +1 (509) 783-9850, Option 4
Web: www.bruker.com/hhxf
Address: 415 N. Quay Street
 Kennewick, WA 99336 USA

1.1 Installation

- Hardware requirements** To use Artax, a PC with Windows 7, 8.1 or 10 operating system and a USB port is required.
- Instructions** *To install Artax –*

Step	Action
1	Insert the USB flash drive.
2	Navigate to the USB drive (part number 160.0211).
3	Double-click Artax – 8-4xx.exe and follow the directions on the screen.

1.2 Program Start



- Description** This section describes how to start Artax and connect to an instrument.
- Note:** A default printer must be accessible or Artax may not start correctly.
- Logging in** *To log into Artax –*

Step	Action	Result
1	On the computer, start Artax by double-clicking the icon.	The Password dialog box is displayed.
2	Enter a User name (test) and Password (test), and press Enter or click Ok .	The Artax window is displayed and OFF-LINE is displayed in the status bar at the bottom.

Connecting to an instrument

To connect to an instrument –

Step	Action	Result
1	Power on the PC and the instrument.	
2	Connect the instrument to the PC with a USB cable (Bruker part number 160.0012) or through Wi-Fi. To enable USB or set up the Wi-Fi, see the <i>Supervisor Manual</i> , Bruker document 030.0113.	
3	From the Artax window, click Device (left side of menu bar).	The dropdown menu is displayed.
4	Click Connect .	The Connect dialog box is displayed.
5	Click the instrument name and Connect .	For the TRACER 5, the Information popup with instrument information is displayed.
6	For the TRACER 5, click OK .	Both dialog boxes are removed. Below the Artax toolbar, as data is retrieved from the instrument, a colored bar is displayed: first yellow, then orange, and finally green. ON-LINE is displayed in the status bar at the bottom. Artax logs into the instrument and the Ready to Test screen is displayed there.



Failure to connect

If Artax is unable to connect to the instrument, either **Error transmission to device** or **Time out** is displayed. Click **OK**.

Artax cannot connect with an instrument if –


- The instrument is not powered on and logged into.
- The USB connection cable between the computer and instrument is not plugged in.
- The computer hardware or the instrument has a malfunction.

2 User Interface

Description This section describes the organization of Artax options.

Window structure The Artax program window contains the –

- Menu bar.
- Toolbar.
- Active panels.
- Status bar.




2.1 Menu Structure

Description The Artax menu bar contains **File, Device, Measurement, Analyze, Spectrum, Project, Options, Export, User, and Help**. Their submenus are described below.


2.1.1 File Menu

Description This section describes options found under **File** in the menu bar. Note that during a measurement these options are not available.



Option	Description
Open Project, Open ROI, Open Spectrum	Opens a dialog box through which a project, ROI (region of interest), or spectrum/spectra can be selected and opened. The most recent directory and files are listed.
Reopen	Displays a list of the last modified spectra or project files. Click to open the file. The submenu option Clear deletes all listed items.
Save Project As, Save ROI As, Save Spectrum As	Opens a dialog box through which a file can be saved. The most recent directory and files are listed.
Exit	Closes Artax. If an assay is in process, the assay is stopped. The message Project was changed! Exit Program? is displayed if either of the following are unsaved – <ul style="list-style-type: none"> • A modified open project. • Measurement spectra in the Spectrum panel.


2.1.2 Device Menu



Description This section describes options found under **Device** in the menu bar.

Option	Description
Connect	Once the PC is physically connected to an instrument via USB cable or Wi-Fi, displays the Connect dialog box through which Artax connects via software. See Program Start on page 5.
Disconnect	Disconnects Artax from an instrument.


2.1.3 Measurement Menu



Description This section describes options found under **Measurement** in the menu bar.

Option	Description
Start	If the PC is connected to the instrument, starts a measurement.
Stop	During a measurement, stops a measurement.
Method	Opens the Method Editor dialog box through which methods can be selected, added, replaced, and removed. See Working with Methods on page 39.


2.1.4 Analyze Menu



Description This section describes options found under **Analyze** in the menu bar.

Option	Description
Accumulate Spectra	Sums – <ul style="list-style-type: none"> Channels using all spectra of one selected node within a project data file. Real and live measuring times. A new node, Accu_Points or Accu_Spectra , is added to the project.
Evaluate Results	Calculates the net intensities via deconvolution for all spectra within one selected node of a project file.
Evaluation	Calculates net intensities using selected deconvolution method.
Match	Opens the Match dialog box through which a spectrum can be compared with saved spectra. See Comparing Spectra on page 47.
Periodic Table	If a spectrum is open, displays the Periodic Table of the Elements dialog box. See Qualitative Analysis on page 49.

2.1.5 Spectrum Menu




Description

This section describes options found under **Spectrum** in the menu bar. The functions of the **Spectrum** menu, except **Delete All ROI**, are not available during measurements.



The *active spectrum* is the currently selected spectrum that can be analyzed. See **Working with Spectra** on page 23.


When measurements are taken, spectra are added, or spectra are subtracted, newly created files are given the temporary names of Measured_n, Add_n, and Sub_n, respectively.

Option	Description
C = A + B	Displays Select Spectra dialog box to indicate two spectra to add. See Adding and Subtracting Spectra on page 45.
C = A - B	Displays Select Spectra dialog box to indicate a spectrum to subtract from another. See Addition and Subtraction of Spectra on page 45.
A = F * A	Displays Factor for multiplication dialog box to enter a multiplication factor to scale the active spectrum. See Scaling a Spectrum on page 46.
Normalize	Scales all displayed spectra to the active spectrum to the channel at the cursor location.
Smooth	Smooths the active spectrum via a polynomial filter of second order. Smoothed values replace original data.
Counts	Opens the Counts dialog box. The cursor can be moved one channel at a time and the counts in a channel can be changed. It also displays the energy and counts of the current channel.




Description,
cont.

Option	Description
Degrade	Opens the FWHM Mn K Alpha dialog box to increases the full width at half maximum (FWHM) of each line of a spectrum. <div data-bbox="794 411 1362 579" style="border: 1px solid black; padding: 5px; margin: 10px 0;">  </div>
Degrade All Reference	Opens the Degrade All Reference dialog box through which the source and destination folder can be entered. Specify resolution (FWHM). Degrades all .spx spectra to the resolution specified and stores them in the destination folder. <div data-bbox="818 737 1338 1010" style="border: 1px solid black; padding: 5px; margin: 10px 0;">  </div>
Close Spectrum	Closes the active spectrum.
Close All Spectra	Closes all open spectra.
Delete ROI	Removes the selected ROI.
Delete All ROI	Removes all ROIs.
Change Parameter	Not currently supported and may give incorrect results.



2.1.6 Project Menu




Description

This section describes options found under **Project** in the menu bar when the project nodes panel on the left of the window is right clicked.

Option	Description
New Project	Creates a new project.
Close Project	Closes a project. Note: this option does not save changes.
Display Spectra	Displays all spectra of the selected node, unless Options -> Only one Spectrum is checked.
Add Spectra	Adds all displayed spectra to the project.
Add Picture	Not currently available.
Add Clipboard	Adds objects on the clipboard to the project.
Add Node	Adds to the project a node.
Remove Object	Removes the selected object from the project.

2.1.7 Options Menu




Description

This section describes options found under **Options** in the menu bar and when the Spectrum panel is right clicked.

Option	Description when checked
Filled Spectrum	Displays a colored area under the line plot.
Filled Background	Displays a colored area under the background radiation curve. Do not set both the spectrum and background to Filled . The area plot of the filled spectrum would occlude the area plot of the background.
Counts per Second	Displays the Y-axis scale's units in counts per second.
Only one Spectrum	Displays only one spectrum. If a new measurement is started or if a new spectrum data file is opened, the Spectrum panel is cleared.
Channel Scale	Displays the X-axis scale's units in channels.
Acoustic signal when finished	Not currently available.
Logarithmic Scale	Displays Y-axis values as a logarithmic scale.
Automatic Scale	During a measurement, adjusts the scale to the maximum amplitude in the active spectrum, excluding the zero peak.
Corrected data curve	Displays a data curve corrected for escape peaks and sum peaks. Corrects only a data curve placed within a project.
Background data curve	Displays the background data curve.
Difference curve	Displays the difference between the original spectrum and the deconvoluted spectrum.
Short Result List	Displays the shortened intensity list of deconvolution results to the right of the spectrum.
Display Properties	Opens the Display Properties dialog box which defines how spectra are displayed. See Display Options on page 25.
Report Designer	Not currently available.
Report Options	Not currently available.

2.1.8 Export Menu




Description

This section describes options found under **Export** in the menu bar.


Option	Description										
Spectrum	Exports spectrum data as a text file (.txt).										
All Spectra	Exports all spectra contained in a node as a text file.										
Result	Exports deconvolution results as a text file.										
All Results	Exports all results contained in a node as a text file.										
Results to Excel	For each analyzed element, exports the results of a line-scan or mapping to an .xlsx file that can be opened by Microsoft Excel.										
All Picture	Not currently available.										
Options	Opens the Export Options dialog box which controls what is included in the text file created from Export options Spectrum through Results to Excel .										
Copy Graphic	<table border="1"> <thead> <tr> <th>Active Tab</th> <th>Exports to the Windows Clipboard</th> </tr> </thead> <tbody> <tr> <td>Spectrum</td> <td>Spectra displayed in the Spectrum panel.</td> </tr> <tr> <td>Mapping Results</td> <td>Not currently available.</td> </tr> </tbody> </table>	Active Tab	Exports to the Windows Clipboard	Spectrum	Spectra displayed in the Spectrum panel.	Mapping Results	Not currently available.				
Active Tab	Exports to the Windows Clipboard										
Spectrum	Spectra displayed in the Spectrum panel.										
Mapping Results	Not currently available.										
Copy Result	Copy results in table format. <table border="1"> <thead> <tr> <th>Active Tab</th> <th>Exports to the Windows Clipboard</th> </tr> </thead> <tbody> <tr> <td>Spectrum</td> <td>Summary of the deconvolution results.</td> </tr> <tr> <td>Results</td> <td>Short result list.</td> </tr> <tr> <td>Mapping Results</td> <td>Not currently available.</td> </tr> <tr> <td>Match Results</td> <td>Results table.</td> </tr> </tbody> </table>	Active Tab	Exports to the Windows Clipboard	Spectrum	Summary of the deconvolution results.	Results	Short result list.	Mapping Results	Not currently available.	Match Results	Results table.
Active Tab	Exports to the Windows Clipboard										
Spectrum	Summary of the deconvolution results.										
Results	Short result list.										
Mapping Results	Not currently available.										
Match Results	Results table.										
Copy ROI	Exports the ROI calculation results to the Windows clipboard.										
Copy Picture	Not currently available.										

2.1.9 User Menu



Description Options under **User** are not currently applicable.

2.1.10 Help Menu














Description This section describes options found under **Help** in the menu bar.

Option	Description
Help	Currently, this <i>Artax User Manual</i> contains the most updated material. Do not use the Help option.
About	Displays Artax version and copyright information. To remove it, press Esc on the keyboard or click the message.



2.2 Toolbar Structure

Description The toolbar displays icons to easily access commonly used options.

Icons Note that icons, except for the stop icon, are not available during a measurement.

Icon	Description									
	Opens the Open Project dialog box through which a project can be selected and opened. The most recent directory and files are listed.									
	Opens the Save Project dialog box through which a project can be saved. The most recent directory and files are listed.									
	Displays the name of the active spectrum. To activate a different spectrum from currently open spectra, click the down pointing caret and click a spectrum name. That spectrum is now active and its name is displayed in the field.									
	Removes all spectra from the Spectrum panel.									
	Switches the Y-axis between linear and logarithmic scale.									
	<p>Activates zoom mode. <i>To zoom in on part of a spectrum –</i></p> <table border="1" data-bbox="716 1073 1516 1367"> <thead> <tr> <th data-bbox="716 1073 797 1106">Step</th> <th data-bbox="797 1073 1122 1106">Action</th> <th data-bbox="1122 1073 1516 1106">Result</th> </tr> </thead> <tbody> <tr> <td data-bbox="716 1106 797 1220">1</td> <td data-bbox="797 1106 1122 1220">Click the zoom icon and move the mouse cursor to the Spectrum panel.</td> <td data-bbox="1122 1106 1516 1220">The mouse cursor displays as a magnifying glass.</td> </tr> <tr> <td data-bbox="716 1220 797 1367">2</td> <td data-bbox="797 1220 1122 1367">Left click, drag to define a rectangle, and release.</td> <td data-bbox="1122 1220 1516 1367">The area within the rectangle is enlarged to fill the Spectrum panel and the mouse cursor reverts to an arrow.</td> </tr> </tbody> </table>	Step	Action	Result	1	Click the zoom icon and move the mouse cursor to the Spectrum panel.	The mouse cursor displays as a magnifying glass.	2	Left click, drag to define a rectangle, and release.	The area within the rectangle is enlarged to fill the Spectrum panel and the mouse cursor reverts to an arrow.
Step	Action	Result								
1	Click the zoom icon and move the mouse cursor to the Spectrum panel.	The mouse cursor displays as a magnifying glass.								
2	Left click, drag to define a rectangle, and release.	The area within the rectangle is enlarged to fill the Spectrum panel and the mouse cursor reverts to an arrow.								
	If the PC is connected to the instrument, starts a measurement.									
	Stops a measurement.									
	Determines the length in seconds of a measurement. To change the duration, click the up pointing or down pointing triangles.									
	Opens the Method Editor dialog box through which methods can be selected, added, replaced, or removed. See Working with Methods on page 39.									
	Lists available methods. A method name may be typed in the field.									

Icons, cont.

Icon	Description
	Displays the Periodic Table of the Elements dialog box. See Qualitative Analysis on page 49.
	Deconvolutes the active spectrum and calculates element intensities.

Radiation icons

During a measurement, two radiation warning icons are displayed in the toolbar –



3 Working with Projects

Description Spectra and analysis results are saved in project files and are automatically created when a measurement is performed through Artax. A project file can contain –

- Several spectra.
- Images.
- Analysis result sets.
- Results of deconvolution images.
- Comments.


3.1 Opening a Project

Description An existing project can be opened and analyzed. A project can contain spectra, ROI (region of interest), and calculation results.

Opening a project

To open a project –


Step	Action	Result
1	Click File .	A dropdown list is displayed.
2	Click Open Project . Note: if a project is already open, clicking Open Project will close it without saving changes.	The Open Project dialog box is displayed.



3	Click the down pointing triangle of the Look in field.	A dropdown list of folders is displayed.
4	Navigate the list and click a folder.	The node is displayed in the Name column.
5	Double-click the project name.	The project nodes panel is opened on the left with nodes listed in a tree. Available tabs are displayed above the main panel. The display in the main panel depends on which tab is active. For information on tabs, see Tabs on page 19.

Displaying a spectrum

To display a spectrum, click its name in the node tree and click the **Spectrum** tab.




Editing a node name

To edit a node name, click it once to select it, then slowly click it again. The name changes to an editable field. Note that the names “Project” and “Objects” cannot be changed.

3.2 Tabs


Description

Various tabs, visible under the toolbar and to the right of the project nodes panel, display specific panels.



Tabs

Tab	Displays	When Visible	Reference
Spectrum	Spectra, a shortened results table of the Bayes deconvolution, and an ROI table, depending on options selected under Options .	Always.	Working with Spectra , page 23.
Parameter	Parameters of an individual spectrum.	After starting a measurement or after displaying a spectrum.	Spectra Parameters , page 35.
Results	Complete Bayes-deconvolution table of results.	After executing a deconvolution or after displaying a spectrum containing results.	Deconvolution , page 51.
Project	User - Name of the user under whose login the project was created. Comment - Provides the ability to enter text. To save changes, save the project.	After loading a project file.	Creating a Project , page 22.




Spectrum panel

Parameter panel

No.	Element	Line	Energy/keV	Cycl.	Net	Backgr.	Sigma	Chi	Conc.	SigmaC	LLD
1	K	K12	3.314	9	1534	3993	98	0.80	0.000	0.000	0.000
2	Ca	K12	3.692	9	1030	2878	82	0.56	0.000	0.000	0.000
3	Cr	K12	5.415	9	237	1037	48	1.87	0.000	0.000	0.000
4	Mn	K12	5.900	9	1129	1066	57	1.68	0.000	0.000	0.000
5	Fe	K12	6.405	9	24559	1056	163	4.89	0.000	0.000	0.000
6	Cu	K12	8.046	9	1316	828	55	1.44	0.000	0.000	0.000
7	Zn	K12	8.637	9	863	943	52	0.84	0.000	0.000	0.000
8	Ga	K12	9.251	9	375	912	47	1.00	0.000	0.000	0.000
9	Ga	L1	1.098	9	-1	6745	116	2.11	0.000	0.000	0.000
10	As	K12	10.543	9	545	891	48	0.85	0.000	0.000	0.000
11	As	L1	1.282	9	1	6260	112	0.68	0.000	0.000	0.000
12	Kr	K12	12.648	9	599	1251	56	2.91	0.000	0.000	0.000
13	Kr	L1	1.585	9	40	6000	110	2.24	0.000	0.000	0.000
14	Rb	K12	13.396	9	7708	1610	105	2.15	0.000	0.000	0.000
15	Rb	L1	1.692	9	169	5930	110	2.00	0.000	0.000	0.000
16	Sr	K12	14.165	9	24653	1600	167	6.94	0.000	0.000	0.000
17	Sr	L1	1.806	9	79	6329	113	0.12	0.000	0.000	0.000
18	Y	K12	14.958	9	2121	1772	75	1.08	0.000	0.000	0.000
19	Y	L1	1.924	9	79	6213	112	0.89	0.000	0.000	0.000
20	Rh	K12	20.216	9	6087	16188	196	2.78	0.000	0.000	0.000
21	Rh	L1	2.697	9	45	5051	101	1.18	0.000	0.000	0.000
22	Pd	K12	21.177	9	54600	22040	314	26.93	0.000	0.000	0.000
23	Pd	L1	2.838	9	616	4779	101	0.98	0.000	0.000	0.000
24	Sn	K12	25.271	9	6041	47655	318	1.73	0.000	0.000	0.000
25	Sn	L1	3.444	9	0	3595	85	0.96	0.000	0.000	0.000
26	Ba	K12	32.194	9	5576	51130	328	1.51	0.000	0.000	0.000
27	Ba	L1	4.466	9	510	1473	59	4.67	0.000	0.000	0.000

Results panel




Project panel

3.3 Creating a Project

Description A project is a collection of measurement spectra. A project can be created before taking measurements and after.

After measurements *To create a project after taking measurements –*

Step	Action	Result
1	Open spectra or take measurements to be included in the new project.	
2	Click Project .	A dropdown list is displayed.
3	Click New Project .	The project nodes panel is labeled New project .
4	In the project nodes panel, click to highlight Objects .	
5	Click Project from the menu bar and Add Node .	The Node dialog box is displayed.
6	Enter a node name and click OK .	The dialog box is removed and the new node is listed under Objects node of the nodes panel.
7	Double click Objects and click a node name.	The name is highlighted.
8	Click Project from the menu bar and Add Spectra .	All open spectra are added to the new project.



Before measurements *To create a project before taking measurements, follow the steps in the preceding table except move step 1 to right after step 3, then resume with step 4.*



3.4 Copying and Saving

Copying *To copy a node to a different object or under another node*, click the node to be copied, drag, and drop it on top of the name of the object or node under which it will go.

Saving To save spectra within projects, click **File -> Save Project As**.

Note: using the **Open Project** or **Exit** option under **File** will not save changes. There is also no message to confirm closing a file without saving.

4 Working with Spectra

Description This section describes types of spectra, and how to display and manipulate them in the Spectrum panel. *To display the Spectrum panel*, select the **Spectrum** tab. Each displayed spectrum has an assigned color.

Note that when a large number of spectra are displayed, the general view worsens and performance slows.

For a description of tabs and their panels, see **Tabs** on page 19.

Spectra types For every measured or imported spectrum displayed, four spectra are managed within Artax –

Spectrum Type	Description
Original data	Data from a measurement of a sample.
Corrected spectrum	The spectrum with the background subtracted.
Calculated background count rate	The background data curve.
Recalculated spectrum from the deconvolution procedure.	See Deconvolution on page 51.


Assay result lines If assay results are saved with the spectrum, then the lines of the analyzed elements are displayed on the Spectrum panel. See **Qualitative Analysis** on page 49.

4.1 Active Spectrum

Active spectrum

The active spectrum is the currently selected spectrum that can be analyzed.

- It is brought to the top of all other displayed spectra.
- Its name is displayed in the toolbar field.
- Its color is indicated in a box to the left of the name in the toolbar field.
- Its name in the project nodes panel is highlighted.



Selecting a spectrum

To make a spectrum active –


Step	Action	Result
1	Click the down pointing triangle on the left field in the toolbar.	A dropdown list is displayed.
2	Click a spectrum name.	The selected spectrum is active.

OR click the spectrum name in the project nodes panel.

4.2 Display Options

Description What is displayed in the Spectrum panel depends on options selected under **Options -> Display Properties**.

Spectra frame Options in the **Spectra** frame affect how spectra are displayed.



Option	Description
Sequence of colors	As each spectrum opens it is assigned a color according to the sequence defined here. To change a color box , click it. The Color dialog box is displayed. Click a color and OK . The selected color box reflects the choice. To change the order of colors , hold the Alt key, click and drag a color box to a new location, and release. The two boxes are exchanged.
Interpolate Type	None – Displays the spectrum with constant count rate across the width of the channel. Linear – The spectrum curve is linearly interpolated between channels.
Background box	Color of the background curve.
Reconstructed spectrum box	Color of the reconstructed spectrum after deconvolution.

Customizing colors

The **Color** dialog box provides the option to create custom colors. *To create a custom color –*

Step	Action	Result
1	From the Display Properties dialog box, click a color box.	The Color dialog box is displayed.
2	Click Define Custom Colors .	The dialog box expands.
3	Select a color by clicking in the color gradient or editing the values in the numeric fields.	The Color box and numeric fields reflect the color choice.
4	On the right side of the dialog box, slide the left pointing triangle up or down as needed.	The Color box and numeric fields reflect the color choice.
5	Click Add to Custom Colors .	The custom color is displayed with the Custom colors on the lower left side of the dialog box.
6	Click the custom color to apply to the previously-selected Display Properties option and click OK .	The Color dialog box is removed and the previously-selected Display Properties option color box reflects the change.







Chart frame Options in the **Chart** frame affect what is displayed in the Spectrum panel.




Option	Description
Color box	Displays the Color dialog box. To change the chart background color , click the color box, click a color on the displayed Color dialog box and click OK , and click OK on the Display Properties dialog box.
Axes Font	Displays the Font dialog box. To alter the text format of axes markers , click desired options and OK .
Marker Font	See Axes Font above.
Abscissa (low) This option, with the following two, define spectrum panel boundaries; how much of the spectrum is displayed.	Displays the field value as the <i>smallest</i> value of the <i>energy</i> scale. <ol style="list-style-type: none"> 1. Open a spectrum. 2. Open Display Properties and check the Abscissa (low) box. 3. Change the value and click OK. 4. Execute the change by double clicking the spectrum panel.
Abscissa (high)	Displays the field value as the <i>highest</i> value of the <i>energy</i> scale. See the steps for Abscissa (low) .
Ordinate (high)	Displays the field value as the <i>highest</i> value of the <i>impulse</i> scale. See the steps for Abscissa (low) .
Series line width	Defines the line weight of all displayed spectra.

Grid frame Options in the **Grid** frame affect how the grid of the spectra is displayed.



Option	Description
None	No grid is visible.
Long ticks	Tick marks are visible where axes are labeled.
Short ticks	Tick marks are visible at axis marks.
Style	Displays a dropdown list of line style options: points, dot lines, horizontal dot lines, horizontal lines, lines, points, vertical dot lines, and vertical lines.
Color box	Displays the Color dialog box. To change the grid color , click a color and OK .

Cursor frame Options in the **Cursor** frame affect how the cursor is displayed.



Option	Description
Vertical	Vertical bar cursor.
Both	Vertical bar and horizontal bar cursor.
None	No cursor.
Follow	When Both is selected and the cursor is moved, the horizontal bar follows the active spectrum curve.
Free	When Both is selected and the cursor is moved, the horizontal bar moves independently of the spectrum.
Color box	Displays the Color dialog box. <i>To change the cursor color, click a color and OK.</i>

4.3 Manipulating a Spectrum

Description To better view spectrum details, use the following commands.


Mouse and keyboard commands

Command	Mouse	Keyboard
Move cursor	Left button click or drag	Left/right arrow
Compress and stretch vertically	Ctrl + Left button OR Left of the Y-axis, left button	Page up/down or arrow up/down
Compress and stretch horizontally	Ctrl + Left button	Ctrl + left/right arrow
Move horizontally	Below the X-axis, left button	Alt + left/right arrow
Move to the beginning		Home
Move to the end		End
Return to original display	Left button double-click	



Zooming

To zoom in on part of a spectrum –

Step	Action	Result
1	Click the zoom icon,  , and move the mouse cursor to the Spectrum panel.	The mouse cursor displays as a magnifying glass.
2	Left click, drag to define a rectangle, and release.	The area within the rectangle is enlarged to fill the Spectrum panel and the mouse cursor reverts to an arrow.
3	To restore the display , double-click in the Spectrum panel.	The spectrum is restored to its original display.

4.4 Status Bar

Description

The status bar at the bottom of the window displays the values for the point the cursor intersects an active spectrum.

Placing the cursor

To place the cursor on a specific point on a spectrum –

Position the mouse cursor and left click, <i>OR</i> left click and drag the cursor, <i>OR</i> use the left/right arrow keyboard keys.	A vertical line intersecting the spectrum is displayed and data relevant to the intersection point is reflected in the status bar. See below.
---	---

If **Options -> Channel Scale** is selected, the cursor can be placed only onto existing data points.

Status bar fields

The status bar displays the following conditional fields and information –

Field	When Displayed	Description
E	Options -> Channel Scale is <i>off</i> .	Energy in keV.
Ch	Options -> Channel Scale is <i>on</i> .	Channel number (0-2048).
Cnts	Options -> Counts Per Second is <i>off</i> .	Count in the channel.
R	Options -> Counts Per Second is <i>on</i> .	Counts per second (rate).
User	Always.	Identification of logged in user.
READY	The PC is connected to an instrument.	The connection between the PC and the instrument is detected.
OFF-LINE	The PC is not connected to an instrument.	No connection between the PC and the instrument is detected.

5 Regions of Interest (ROIs)

Description An ROI displays, in the ROI panel, the area under a curve between ROI boundaries. Multiple ROIs are defined in a single spectrum but apply to all spectra in the project. With ROIs, portions of the spectrum for the calculation of peak intensities can be defined.

Defining an ROI


To define an ROI –

Step	Action	Result
1	Position the mouse cursor at the left edge of the region to be defined.	
2	a) Press and hold the Shift key, b) left click and hold, c) drag to the opposite edge of the region to be defined, d) and release.	During the operation, the cursor displays as a short vertical line with ROI beneath.

Adjusting an ROI

To adjust left and right edges of an ROI –

Step	Action	Result
1	Press and hold the Shift key and move the mouse cursor over an ROI edge.	The ROI cursor gains left and right pointing arrows.
2	Continue to hold the Shift key, left click and hold the mouse, drag, and release.	The ROI edge moves accordingly.




ROI panel

ROIs

ROI display

How an ROI is displayed depends on options selected under **Options** in the menu bar, **Display Properties, ROI** frame.



Option	Description
Foreground	Not currently available.
Background	ROIs are displayed as shaded areas, each labeled with a name. The shade color is determined by the Color box.
Color box	Displays the Color dialog box. To change the color for new ROIs , click a color and OK . For details, see Display Options on page 25. This does not change the color of existing ROIs.

ROI panel

At the bottom of the window, the ROI panel is displayed. It contains data for all ROIs in the active spectrum. When a different spectrum is made active, the ROIs stay the same, but the data in the ROI panel are recalculated.

Field	Description
No.	Consecutive number based on creation order.
Start/keV	Beginning of the ROI.
End/keV	End of the ROI.
Name	If labeled through the periodic table, a number or element line name.
Gross	Calculated count rate.
Net	Gross minus background count rate.
Backgr.	Background count rate.
Spectrum	Names of spectra containing ROIs.

Copying

To copy data from the ROI panel to the Windows clipboard –

From the menu bar, click **Export -> Copy ROI**.

OR

Right click in the ROI panel and left-click **Copy ROI**.



Naming ROIs are numbered consecutively, 1 to N. If elements are labeled on the spectrum, then the ROI is assigned a standard name of the element symbol and line type (e.g., Cu_K).

Renaming *To rename an ROI –*

Step	Action	Result
1	In the ROI panel, click an ROI.	The row is highlighted.
2	Click Name .	The ROI Name dialog box is displayed.
3	Enter a new name and press Enter or click OK .	The dialog box is removed and the name in the ROI panel reflects the change.

Deleting *To delete a single ROI using the Spectrum panel –*

Step	Action	Result
1	In the Spectrum panel, double-click inside an ROI. (If spectra are expanded, first double-click to resize the display.)	The ROI is bordered by dashed lines.
2	Press Delete on the keyboard.	The ROI is removed from the Spectrum panel and the ROI panel.

To delete a single ROI using the ROI panel –

Step	Action	Result
1	In the ROI panel, click an ROI.	The ROI is highlighted.
2	From the menu bar, click Spectrum . <i>OR</i> Right click the ROI in the ROI panel.	A dropdown list is displayed.
3	Click Delete ROI .	The ROI is removed from the Spectrum panel and the ROI panel.

To delete all ROIs –

Step	Action	Result
1	From the menu bar, click Spectrum . <i>OR</i> Right click in the ROI panel.	A dropdown list is displayed.
2	Click Delete All ROI .	All ROIs are removed from the Spectrum panel and the ROI panel.

Saving

To save an ROI –

Step	Action	Result
1	From the menu bar, click File .	A dropdown list is displayed.
2	Click Save ROI As .	The Save ROI List dialog box is displayed.
3	To select a folder , click in the Save in field.	A dropdown list is displayed.
4	Click the desired folder.	Folders and files contained in the selected folder are displayed.
5	Enter a name in the File name field and click Save or press Enter.	The ROIs are saved in a file and the dialog box is removed.

Opening

To open a saved ROI file –

Step	Action	Result
1	From the menu bar, click File .	A dropdown list is displayed.
2	Click Open ROI .	The Open ROI List dialog box is displayed.
3	In the Look in field, click the down pointing triangle to navigate to the desired folder.	The desired file is listed in the Name field.
4	Double-click the file name. <i>OR</i> Click the file name and click Open or press the Enter key.	ROIs from the file are added to the Spectrum panel and the ROI panel.



6 Spectra Parameters

Description The Parameter panel, visible when the **Parameter** tab is selected, displays measurement parameters for the active spectrum.

Spectrum parameters

Field	Description
Date/Time	Date and time the measurement started. Format depends on operating system settings.
Channel count	Number of channels in the spectrum.
Energy linear	The linear parameter A of the Energy-to-channel equation ($E=A*x+B$).
Fano factor	Describes the fluctuation of the primary ionization at the detector.
Counts Per Sec	Counts per second over the entire spectrum.

Excitation parameters

Field	Description
Anode	Anode material of the X-ray tube.
Filter	Filter material.
Atmosphere	Milieu of the measuring chamber: air, He, or vacuum.
Pressure	Pressure in the measuring chamber in mbar.
Temp	Temperature in Celsius of the measuring chamber.

Measurement parameters

Field	Description
Input counts/sec	All counts, including valid, double, escape, and more.
Output counts/sec	Only valid counts.
Dead time	Percent of time counting rejected invalid counts.

X-ray Generator parameters

Field	Description
Voltage	High voltage of X-ray tube.
Current	X-ray tube current.


Method parameters

These values are set in the **Method Editor**. See **Working with Methods** on page 39.

7 Exporting

Description Measurement data can be exported for further analysis to standard programs (e.g., Microsoft Excel) and other evaluation software.

Settings *To define what data is exported*, check settings in the **Export Options** dialog box, under **Export -> Options**.



Spectra frame

Option	Description
Only counts	Only recorded pulse counts.
Channels and counts	Channel numbers and recorded pulse counts.
Energy and counts	Energy values and recorded pulse counts.
Number per Line	Number of measurements per line of text.
With header	Measurements are preceded by the eV/channel and zero offset. Exported spectra files (.txt files) can be reimported, if With Header was checked when they were exported. <i>To import</i> , use File -> Open Spectrum....

Results frame

Option	Description
Net area	Calculated net intensities.
Concentration	Not currently available.

Active spectrum data

To export original measurements of the active spectrum –

Step	Action	Result
1	From the menu bar, click Export -> Spectrum .	The Save Spectrum as ANSI File dialog box is displayed.
2	Navigate the Save in field to locate the folder to export to.	The folder name is displayed in the Save in field.
3	Provide a name in the File name field and click Save .	Spectrum data, as specified in Export Options dialog box, are saved in a file in the selected folder.

All spectra
data

To export all spectra from a project node as text files –

Step	Action	Result
1	Select a spectrum or node from the project nodes panel.	Spectra are displayed.
2	Click Export -> All Spectra .	The Save Spectra from Node as ANSI File dialog box is displayed.
3	Navigate the Save in field to locate the folder to export to. The original file path is recommended.	The folder name is displayed in the Save in field.
4	Provide a name in the File name field and click Save . That name is used as a prefix for each created file. The original file name with extension .txt is recommended.	Each spectrum's data are saved in an individual file in the selected folder.

Current results

To export results of the deconvolution and intensities (counts) of the active spectrum –

Step	Action	Result
1	Click Export -> Result .	The Save Results as ANSI File dialog box is displayed.
2	Navigate the Save in field to locate the folder to export to.	The folder name is displayed in the Save in field.
3	Provide a name in the File name field and click Save .	Data, as specified in Export Options dialog box, are saved in a file in the selected folder.

All results

To export all results from a project file node to a text file in .csv format –

Step	Action	Result
1	Select a spectrum or node from the project nodes panel.	Spectra are displayed.
2	Click Export -> All Results .	The Save Results from Node as ANSI File dialog box is displayed.
3	Navigate the Save in field to locate the folder to export to. The original file path is recommended.	The folder name is displayed in the Save in field.
4	Provide a name in the File name field. The original file name with extension .csv is recommended.	Data, as specified in Export Options dialog box, are saved in a file in the selected folder.

To Excel

To export all results from a project file node to an Excel file –


Step	Action	Result
1	Select a spectrum or node from the project nodes panel.	Spectra are displayed.
2	Identify element lines in the spectra using Auto Ident in the Periodic Table of the Elements . See Qualitative Analysis on page 49.	Element peaks are identified.
3	Click Analyze -> Evaluate Results . The more spectra to analyze, the longer the analysis.	The intensities of the selected lines for all spectra are determined and stored in memory. A progress bar may be displayed.
4	Click Export -> Results to Excel .	The Save As dialog box is displayed.
5	Provide a file name with an .xlsx extension and click Save .	The dialog box is removed and the results are saved as an Excel file.


8 Working with Methods

Description Parameters for measuring and evaluating spectra data are specified within one method. An unrestricted number of methods, classified by name, can be saved.

Selecting a method


To select a different method –


Step	Action	Result
1	Click the down pointing triangle to the right of the Edit methods icon ()	A dropdown list of method names is displayed.
2	Click a method name.	The list is removed, the selected method name is displayed in the field, and relevant parameters are updated.



8.1 Method Editor

Description The **Method Editor** is a dialog box for creating new methods and editing existing ones.

To display the Method Editor, click the **Edit methods** icon () from the toolbar OR, from the menu bar, click **Measurement -> Method**. **To access different panels**, click the appropriate tab.




General fields

Option	Description
Comment	<p>Comments entered for the currently selected method.</p> <p>When the method name is displayed in the toolbar method box and the mouse cursor hovers over the name, text previously entered in this Comment field is displayed.</p>
Name	Currently selected method.
Add	<p>To add a new method, enter a unique name in the Name field and click Add. The new name is displayed in the list box.</p>
Replace	<p>To change parameters of the selected method, click the method name, edit field values, and click Replace.</p>
Remove	Deletes the selected method.
Automatic evaluation	<p>When a measurement completes, deconvolution (blue curve on the Spectrum panel) and intensities (data displayed on the Results panel) are calculated. Detected elements are labeled and saved with the spectrum. The Bayesian deconvolution curve automatically changes when a new element is added or subtracted via the periodic table. See Qualitative Analysis on page 49.</p> <p>This option is available if the following are selected –</p> <ul style="list-style-type: none"> • Under the Deconvolution tab, Bayes OR Profile bayes (normal Fit) AND • Under the Identification tab, Line markers OR Automatic. <p>Note: Intensities results are displayed in the Results panel only if the spectrum is in a project.</p>
Data and time	The date and time of the last method modification.
Ok	Stores all changes in the database.

8.2 Measurement Panel

Description




Option	Description
Voltage	X-ray tube voltage.
Current	X-ray tube current.
Max. Current	When clicked, displays in the Current field the maximum current possible at the voltage setting.
Filter	A dropdown list of possible material filters.
Time	Duration of a measurement.
Select Illumination	From the list of available illuminations, select one as a starting point.
Auto-save PDZ file(s) in location	When this box is checked, saves the .pdz file in the folder or file specified in the boxes below.
Folder	Folder in which the PDZ file is saved.
File Name	PDZ file name with _n.pdz automatically appended, where n is a unique number. The default name is Measured.

TRACER 5 options

When connected to a TRACER 5, the following additional options are available.





Option	Description
Atmosphere	Atmospheric measurement environments: Air – Ambient atmosphere. Vacuum – setting for use with a vacuum pump to eliminate air between the sample and detector. Flush – setting for use with a Helium, or other gas, flush to replace air between the sample and detector.
Collimator	Collimators that can be manually installed.
Manual filter	Filters that can be manually installed.

8.3 Corrections Panel




Option	Description
Escape	Corrects the spectrum for escape peaks.
Shelf	Corrects the spectrum for background signals in the lower energy range (Shelf).
Background	Corrects the background.
Pileup	Ignores sum peaks.
Cycles	For background, tunes background calculations.
Start	For background, lowest energy boundary for correction.
End	For background, highest energy boundary for correction.

8.4 Identification Panel


Method Name	Description
Line markers	<p>Measured elements within the spectrum must be identified via the Periodic Table of the Elements dialog box, via the toolbar icon () or Analyze -> Periodic Table. See Qualitative Analysis on page 49.</p> 
Preset list	<p>This option is recommended when the collected spectra are known to contain certain elements. To select required elements, click element buttons on the periodic table.</p> <p>To add elements to the periodic table from the active spectrum, if the elements there are labeled, click Get Elements.</p> 
Automatic	<p>Identifies elements within the spectrum. To exclude elements known to not be present, click element buttons on the periodic table.</p>  <p>To exclude elements impossible to measure or rare, click Standard.</p>

8.5 Deconvolution Panel


Option	Description
Bayes	Performs Bayesian deconvolution of spectrum. <div style="border: 1px solid black; padding: 5px; margin: 10px 0;">  </div>
Profile bayes (normal fit)	Not currently available.
Profile bayes (optimized fit)	Not currently available.

8.6 Qualification Panel

Not currently available.



8.7 PDZ Options Panel



The PDZ version can be 24 or 25.

9 Spectra Processing

Description Spectra can be added, subtracted, scaled, compared, and saved.

Saving New spectra can be saved as part of a project or as an individual spectrum. Be sure to use a new file name so another file is not overwritten.

*To add a new spectrum to the current project, click **Project -> Add Spectra**.*

*To save a new spectrum as an individual spectrum rather than part of a project, click **File -> Save Spectrum As**.*



9.1 Adding and Subtracting Spectra

Description Two spectra can be added, or one spectrum subtracted from another.

Adding or subtracting

To add two spectra or subtract one spectrum from another –

Step	Action	Result
1	Open two or more spectra.	
2	From the menu bar, click Spectrum -> C = A + B OR C = A - B .	The Select Spectra dialog box is displayed.
3	Under Spectrum A , click the down pointing caret to display a list of open spectra and click a spectrum name. Parameter values from spectrum A will apply to the new spectrum.	The list is removed and the selected spectrum name is displayed in the box.
4	Under Spectrum A , select the appropriate option to be added or subtracted. Original values – The original spectrum. Corrected values – After the background is subtracted. Background – Just the background. Fitted values – The deconvoluted spectrum.	
5	Under Spectrum B , click the down pointing caret to display a list of open spectra and click a spectrum name.	The list is removed and the selected spectrum name is displayed in the box.
6	Under Spectrum B , select an appropriate option and click Ok .	The dialog box is removed; the Spectrum panel displays the new spectrum, Add_n or Sub_n (n being a consecutive number); and the toolbar indicates that it is the active spectrum.




9.2 Scaling a Spectrum

Description Scaling is multiplying the pulse counts, which affects the vertical scale.

Scaling *To scale the active spectrum –*

Step	Action	Result
1	Click Spectrum -> A = F * A .	The Factor for multiplication dialog box is displayed.
2	Type a scaling factor and click OK .	All pulse counts and measurement times of the active spectrum are multiplied by the entered factor.



Normalizing This type of scaling –

- At the spectrum cursor location, gives all open spectra the same value as the active spectrum.
- Recalculates all open spectra by the scaling factor determined by the relation of the pulse counts at the spectrum cursor.


To normalize all spectra in the Spectrum panel to the active spectrum, place the spectrum cursor at the desired location and click **Spectrum** -> **Normalize**. Open spectra are scaled accordingly.


9.3 Comparing Spectra

Description An unknown spectrum can be identified by making it the active spectrum and comparing it to known spectra in project files. The comparison calculates matches in percent.

Searching for matches *To search for matches with the active spectrum –*

Step	Action	Result
1	From the menu bar, click Analyze -> Match .	The Match dialog box is displayed.
2	<i>To select a directory in which to search for matching spectra</i> , click the folder icon to the right of the Search in field.	The Browse For Folder dialog box is displayed.
3	Navigate and select the desired folder, then click OK .	The directory name is displayed in the Search in field.
4	In the Start energy field, specify the lower boundary of the energy range in which to compare the spectra.	The field reflects the change.
5	In the End Energy field, specify the upper boundary of the energy range in which to compare the spectra.	The field reflects the change.
6	In the Min. correlation field, specify the minimum percent of matching concentrations for spectra to be a match.	The field reflects the change.
7	In the Number of hits field, specify the maximum number of matches displayed.	The field reflects the change.
8	Click Start .	The status bar of the dialog box displays the number of files found, the number spectra being compared, and the count of matches found so far. When the comparisons end, the dialog box is removed. If matches were found, the Match Results panel is displayed.





Match Results panel

The Match Results panel displays the following columns –

Column Name	Description
Correlation/%	Closeness of match in percent.
Spectrum	Spectrum name.
Date	Date and time the spectrum was recorded.
Filename	Directory and filename of the project file or of the individual spectrum file.


Active match


To select a spectrum match to be active, double click its row in the Match Results panel. The Spectrum panel is displayed and the selected spectrum is active.

Removing the panel

The Match Results panel remains available until –
 A comparison is executed that results in zero hits.
OR
 Artax is exited.

10 Qualitative Analysis


Description Elements to be charted in the Spectrum panel can be specified from the **Periodic Table of the Elements**, which is a dialog box accessible by clicking –
 The toolbar icon ()
 OR
 The menu bar option **Analyze -> Periodic Table**.



10.1 Display Options

Element name and number Element buttons are displayed in periodic table format. **To display an element button's element name and atomic number**, hover the mouse cursor over it.

Element details **To display details of a specific element**, right click the element button. The **Atomic Data** box is displayed, which contains atomic mass, density, line series energies, and more.



Spectrum display options Depending on boxes checked on the **Periodic Table**, the Spectrum display can include –

- Lines indicating K, L, and/or M series points.
- Flags, with the element symbol and color corresponding to the color of the element button, indicating K, L, and/or M series points.



Option	Description												
K series lines L series lines M series lines	Displays the spectrum K, L, and/or M series lines and/or text for each element highlighted on the Periodic Table , if the Lines box and/or Text box are checked.												
Lines	Displays, with vertical lines, the K, L, and/or M points for each element highlighted on the Periodic Table , if the respective series line boxes are checked.												
Text	Displays, with colored flags indicating the element, the K, L, and/or M points for each element highlighted on the Periodic Table , if the respective series line boxes are checked.												
Color	<p>Changes the color of the –</p> <ul style="list-style-type: none"> • Element button in the Periodic Table. • Element flag and element line on the Spectrum panel. <p>To change the color –</p> <ul style="list-style-type: none"> • Left click and drag a color box from the dialog box and drop it over an element button. <p><i>OR</i></p> <ul style="list-style-type: none"> • Select an element button and press Enter twice repeatedly to scroll through color options. <p>New color assignments remain while the program runs. They cannot be saved.</p>												
New element	<p>To activate a new element –</p> <p>Type the element symbol in this field and press Enter. OR click the up arrow or down arrow to scroll through the elements. OR left click an element button.</p>												
New element arrows	Incrementally selects elements up or down from the selected element button.												
Auto Ident	<p>Depending on checked box options, automatically identifies –</p> <ul style="list-style-type: none"> • In the Periodic Table, elements of the active spectrum. • In the Spectrum panel, element series lines. • In the Spectrum panel, elements symbols associated with series lines. <p>Careful verification of automatic results is recommended.</p> <p>To select elements individually, click an element button. To unselect, click again.</p>												
Clear All	<p>Removes from the –</p> <ul style="list-style-type: none"> • Spectrum panel all element labeling. • Periodic Table all activated element buttons. 												
Window	<p>Lists possible elements according to data at the spectrum cursor location.</p> <table border="1" data-bbox="586 1503 1474 1751"> <thead> <tr> <th data-bbox="586 1503 727 1535">Field</th> <th data-bbox="727 1503 1474 1535">Description</th> </tr> </thead> <tbody> <tr> <td data-bbox="586 1535 727 1608">keV</td> <td data-bbox="727 1535 1474 1608">Adjustable bandwidth for recognizing a line along the X-axis. The default is 0.05 keV.</td> </tr> <tr> <td data-bbox="586 1608 727 1640">OZ</td> <td data-bbox="727 1608 1474 1640">Atomic number.</td> </tr> <tr> <td data-bbox="586 1640 727 1671">Element</td> <td data-bbox="727 1640 1474 1671">Symbol of the elements detected at that point.</td> </tr> <tr> <td data-bbox="586 1671 727 1703">Line</td> <td data-bbox="727 1671 1474 1703">Emission line series type.</td> </tr> <tr> <td data-bbox="586 1703 727 1751">E/keV</td> <td data-bbox="727 1703 1474 1751">Energy in keV.</td> </tr> </tbody> </table> <p>If Profile bayes is chosen as the deconvolution method, the list of possible elements is further restricted based on the included profile files. Element symbols with no corresponding reference data are underlined.</p>	Field	Description	keV	Adjustable bandwidth for recognizing a line along the X-axis. The default is 0.05 keV.	OZ	Atomic number.	Element	Symbol of the elements detected at that point.	Line	Emission line series type.	E/keV	Energy in keV.
Field	Description												
keV	Adjustable bandwidth for recognizing a line along the X-axis. The default is 0.05 keV.												
OZ	Atomic number.												
Element	Symbol of the elements detected at that point.												
Line	Emission line series type.												
E/keV	Energy in keV.												

10.2 Deconvolution


Description From an existing spectrum, deconvolution creates a model spectrum defined by only elements selected on the **Periodic Table of the Elements** dialog box.

The deconvolution process depends on –

- The method selected in the **Method Editor** dialog box, under the **Identification** tab.
- Availability of other spectra results.


Corrected raw counts

To determine corrected raw counts of the active spectrum via deconvolution –

Step	Action	Result
1	Display element labels and series lines on the active spectrum.	
2	Click  on the toolbar or Analyze -> Evaluation on the menu bar.	The deconvolution curve (default blue) is displayed in the Spectrum panel and the method data and results are temporarily stored in memory.
3	To permanently save changes , use File -> Save Spectrum As or add the spectrum to a project file and use File -> Save Project As .	

Differences

To display the difference between original values and deconvolution values, from the menu bar click **Options -> Difference curve**. Results are displayed in a panel below the Spectrum panel.



Appendix A: Registry Database

Description Artax uses the registry database of the operating system intensively to store the position of all windows and dialog boxes and their display parameters. The information is in the registration key HKEY_CURRENT_USER/Software/ Bruker/ARTAX.

Individual's data If several users are registered within the operating system, an individual's data is also stored in the registry. This data is not automatically available after installation of the software or registration of new users.

System parameters and methods Instrument system parameters and available methods are store in the directory **C:\Program Files (x86)\Bruker\Artax\settings\Spectra** with the following naming convention:

Data Type	Storage File Name
Instrument system parameters	program name.ini
Available methods	program name.mth